

THE DELIVERY LINE ADAPTER

The delivery line adapter provides the guidelines for the new bowler to stand on the mat correctly to establish the delivery line for the jack or bowl whether reading from the bank or natural grass.

The metal base plate is placed over the mat 10 centimeters from the leading edge and held firmly in position by forcing the plate downwards for the two stabilizing pins to penetrate into the surface for stabilization of the plate. The topside of the foot pad at the toe and heel areas are lined with Velcro for the stabilizing of the foot pad when placed in position onto the base metal plate. The foot pad itself is also lined with Velcro at the toe and heel areas on both sides of the foot pad incorporating the change over for a right or left hand player by flipping the pad over, the full length of the foot pad is also lined with Velcro which provides stability of the anchored foot. The foot pad can be adjustable for any angle both forehand and backhand.

The centre line if required can be attached to the base plate by velcro as well as the Delivery line in which the Delivery line can be adjusted along the front of the base plate according to the angle required. The diagram below illustrates the layout and attachments.

